

Exo Lab Italia srl Località Boschetto di Pile (AQ), Italy 02294980681 Phone: +393806812563

www.exolabitalia.com

## PRODUCT DETAILS OF CITRUS LIMON PDNVs (ITbio-PDNVs $^{\text{TM}}$ )

## (LYOPHILIZED FORM)

	Citrus limon PDNVs (Plant-Derived NanoVesicles)
	are nano-sized vesicles surrounded by a lipid-enriched
	membrane and containing different natural
	compounds, including vitamins, antioxidants, proteins,
Description	nucleic acids and other metabolites. Exo Lab Italia ®
Description	
	purify and concentrate PDNVs exclusively from Citrus
	limon from Italian organic agriculture. The PDNVs are
	purified and concentrated using ultracentrifugation
	process in a sterile environment.
Origin	Citrus limon
PDNVs Number	2 x 10 <sup>10</sup>
Size	140.3 +/- 1.9nm
pH value	3.14
Form	Lyophilized
Reconstitution	Reconstitute the lyophilized PDNVs with 500-1000 μl
Reconstitution	of sterile H <sub>2</sub> O or PBS
	PDNVs were characterized by concentration and size
Characterization	distribution, pH value, and Tetraspanin-8 (TET-8)
	expression.
Shipping conditions	2-8°C
	ITbio-PDNVs <sup>™</sup> can be stored at +4°C for a short time
	and at -80°C for no more than 1 year. It is
Storage Conditions	recommended to aliquot PDNVs in the desired volume
	and store asap -80°C (active for up to 1 year). Avoid
	repeated freeze-and thaw cycles.
	- For <i>in vitro</i> and <i>in vivo</i> studies
Applications	- Loading of molecules
**	- Drug delivery system
	This product is intended for laboratory research use
Y ( ) 1 1 Y	only. It is not intended for any animal or human
Intended Use	therapeutic use, any human or animal consumption, or
	any diagnostic use.
	Exo Lab Italia srl works in respect of the normative of
BLS 1	BLS 1
	Preparation of reconstituted PDNVs:
	1. Add 500-1000 μl of sterile H <sub>2</sub> O or PBS
	2. Mix by pipetting up and down and by
	vortexing
Handling Procedures	3. Store at 4°C short term (see Storage
	Conditions section)
	,
	4. For long term storage, prepare multiple
	aliquots for use and store at -20°C/-80°C (see



Exo Lab Italia srl Località Boschetto di Pile (AQ), Italy 02294980681

Phone: +393806812563 www.exolabitalia.com

	Storage Conditions section), avoiding freeze/thaw cycles.
	necze/maw cycles.
Material Citations	If use of this material results in a scientific
	publication, please cite the material in the following
	manner:
	ITbio-PDNVs <sup>TM</sup> from <i>Citrus limon</i> by Exo Lab Italia
References	at Tecnopolo d'Abruzzo (AQ) - Italy www.exolabitalia.com
References	Exo Lab Italia ® represents and warrants:
Warranty	(i) any declaration contained in this
	document, particularly in relation to the
	identity of the PDNVs and its
	composition;
	(ii) to identify the applicable legal
	requirements under the legislation on
	protection of workers and environment;
	(iii) to keep track of the quantities of PDNVs
	used and to be able to provide
	documentation (e.g., shipping documents,
	disposal records, process leakage
	information, etc.) that attests the fact that
	records of these quantities are kept;
	(iv) that PDNVs are not supplied by Exo Lab
	Italia ® to the public;
	(v) that the substance will only be handled by internal staff appointed by Exo Lab Italia
	®.
	This product is intended for laboratory research use
	only. It is not intended for any animal or human
	therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is
	prohibited without a license from Exo Lab Italia ®.
	While Exo Lab Italia ® uses reasonable efforts to
	include accurate and up-to-date information on this
	product sheet, Exo Lab Italia ®makes no warranties or
	representations as to its accuracy. Citations from
Disclaimers	scientific literature and patents are provided for informational purposes only. Exo Lab Italia ®does not
Discianners	warrant that such information has been confirmed to
	be accurate or complete and the customer bears the
	sole responsibility of confirming the accuracy and
	completeness of any such information.
	This product is sent on the condition that the customer
	is responsible for and assumes all risk and responsibility in connection with the receipt, handling,
	storage, disposal, and use of the Exo Lab Italia srl
	product including without limitation taking all
	appropriate safety and handling precautions to



Exo Lab Italia srl Località Boschetto di Pile (AQ), Italy 02294980681

Phone: +393806812563 www.exolabitalia.com

minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the Exo Lab Italia ®product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall Exo Lab Italia ®, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, Exo Lab Italia ® is not liable for damages arising from the misidentification or misrepresentation of such materials.

